

---

View Data Log Report**Filter From** : 01-01-2020 00:00 **Filter To** : 09-01-2020 23:59**Records From** : 08-01-2020 16:25 **Records To** : 09-01-2020 15:55**Equipment Name** : HIMEDIA BOD SMH-33 673/663**Temperature (°C)** **Low Alarm** : 0.0 **High Alarm** : 0.0

---

<b>Temperature (°C)</b>	<b><u>Minimum</u></b>	<b><u>Maximum</u></b>	<b><u>Average</u></b>
Ch-1	9.6	27.1	10.9
Ch-2	10.0	27.1	11.3
Ch-3	9.7	27.0	11.2
Ch-4	9.9	27.0	11.6

## View Data Log Report

**Filter From** : 01-01-2020 00:00 **Filter To** : 09-01-2020 23:59

**Records From** : 08-01-2020 16:25 **Records To** : 09-01-2020 15:55

**Equipment Name** : HIMEDIA BOD SMH-33 673/663

**Temperature (°C)** **Low Alarm** : 0.0 **High Alarm** : 0.0

	Date	Time	Set	Temperature (°C)			
Rec. No	dd-MM-yyyy	HH:MM	Temp	Ch-1	Ch-2	Ch-3	Ch-4
1	08-01-2020	16:25	45.0	26.5	26.5	26.6	26.4
2	08-01-2020	16:30	10.0	27.1	27.1	27.0	27.0
3	08-01-2020	16:35	10.0	21.3	22.6	24.0	21.2
4	08-01-2020	16:40	10.0	14.3	16.3	18.7	14.7
5	08-01-2020	16:45	10.0	9.8	11.5	13.7	10.5
6	08-01-2020	16:50	10.0	10.5	11.6	12.6	11.1
7	08-01-2020	16:55	10.0	10.4	11.2	11.7	10.9
8	08-01-2020	17:00	10.0	10.3	11.0	11.2	10.8
9	08-01-2020	17:05	10.0	10.2	10.7	10.9	10.6
10	08-01-2020	17:10	10.0	10.1	10.7	10.7	10.6
11	08-01-2020	17:15	10.0	10.1	10.6	10.6	10.4
12	08-01-2020	17:20	10.0	9.9	10.5	10.4	10.1
13	08-01-2020	17:25	10.0	9.7	10.3	10.4	10.0
14	08-01-2020	17:30	10.0	9.8	10.3	10.2	9.9
15	08-01-2020	17:35	10.0	9.8	10.3	10.2	10.0
16	08-01-2020	17:40	10.0	9.9	10.3	10.3	10.1
17	08-01-2020	17:45	10.0	10.1	10.6	10.4	10.3
18	08-01-2020	17:50	10.0	10.0	10.5	10.4	10.3
19	09-01-2020	10:35	10.0	21.7	21.6	21.6	21.5
20	09-01-2020	10:40	10.0	20.5	21.1	20.8	21.3

## View Data Log Report

Filter From : 01-01-2020 00:00 Filter To : 09-01-2020 23:59

Records From : 08-01-2020 16:25 Records To : 09-01-2020 15:55

Equipment Name : HIMEDIA BOD SMH-33 673/663

Temperature (°C) Low Alarm : 0.0 High Alarm : 0.0

	Date	Time	Set	-----				Temperature (°C)	-----			
Rec. No	dd-MM-yyyy	HH:MM	Temp	Ch-1	Ch-2	Ch-3	Ch-4					
21	09-01-2020	10:45	10.0	15.1	17.1	16.3	18.6					
22	09-01-2020	10:50	10.0	10.6	12.5	11.6	14.8					
23	09-01-2020	10:55	10.0	9.9	11.1	10.3	12.8					
24	09-01-2020	11:00	10.0	9.8	10.6	10.1	11.9					
25	09-01-2020	11:05	10.0	9.8	10.4	10.0	11.4					
26	09-01-2020	11:10	10.0	9.9	10.3	10.0	11.1					
27	09-01-2020	11:15	10.0	10.0	10.2	10.1	11.0					
28	09-01-2020	11:20	10.0	10.2	10.4	10.2	11.0					
29	09-01-2020	11:25	10.0	9.8	10.2	10.0	10.8					
30	09-01-2020	11:30	10.0	10.1	10.3	10.0	10.8					
31	09-01-2020	11:35	10.0	10.0	10.3	10.0	10.7					
32	09-01-2020	11:40	10.0	9.8	10.1	9.8	10.6					
33	09-01-2020	11:45	10.0	9.8	10.1	9.9	10.5					
34	09-01-2020	11:50	10.0	9.9	10.2	9.9	10.6					
35	09-01-2020	11:55	10.0	10.1	10.2	10.1	10.7					
36	09-01-2020	12:00	10.0	9.7	10.2	9.9	10.6					
37	09-01-2020	12:05	10.0	9.7	10.0	9.7	10.5					
38	09-01-2020	12:10	10.0	9.8	10.1	9.9	10.5					
39	09-01-2020	12:15	10.0	10.1	10.2	10.0	10.6					
40	09-01-2020	12:20	10.0	10.0	10.2	10.0	10.5					

## View Data Log Report

Filter From : 01-01-2020 00:00 Filter To : 09-01-2020 23:59

Records From : 08-01-2020 16:25 Records To : 09-01-2020 15:55

Equipment Name : HIMEDIA BOD SMH-33 673/663

Temperature (°C) Low Alarm : 0.0 High Alarm : 0.0

	Date	Time	Set	Temperature (°C)			
				-----			
Rec. No	dd-MM-yyyy	HH:MM	Temp	Ch-1	Ch-2	Ch-3	Ch-4
41	09-01-2020	12:25	10.0	9.6	10.1	9.7	10.4
42	09-01-2020	12:30	10.0	9.8	10.1	9.9	10.5
43	09-01-2020	12:35	10.0	10.0	10.3	10.1	10.6
44	09-01-2020	12:40	10.0	10.1	10.3	10.1	10.6
45	09-01-2020	12:45	10.0	9.8	10.1	9.9	10.5
46	09-01-2020	12:50	10.0	9.8	10.1	9.9	10.5
47	09-01-2020	12:55	10.0	10.0	10.2	9.9	10.5
48	09-01-2020	13:00	10.0	10.2	10.4	10.2	10.6
49	09-01-2020	13:05	10.0	9.8	10.2	10.0	10.5
50	09-01-2020	13:10	10.0	9.7	10.0	9.9	10.6
51	09-01-2020	13:15	10.0	9.9	10.1	9.9	10.6
52	09-01-2020	13:20	10.0	10.0	10.2	10.0	10.5
53	09-01-2020	13:25	10.0	10.1	10.3	10.0	10.6
54	09-01-2020	13:30	10.0	9.8	10.1	9.9	10.5
55	09-01-2020	13:35	10.0	9.7	10.1	9.8	10.5
56	09-01-2020	13:40	10.0	10.5	10.5	10.4	10.9
57	09-01-2020	13:45	10.0	10.0	10.5	10.3	11.0
58	09-01-2020	13:50	10.0	9.9	10.2	10.1	10.7
59	09-01-2020	13:55	10.0	9.8	10.2	9.9	10.6
60	09-01-2020	14:00	10.0	9.6	10.1	9.8	10.6

## View Data Log Report

**Filter From** : 01-01-2020 00:00 **Filter To** : 09-01-2020 23:59

**Records From** : 08-01-2020 16:25 **Records To** : 09-01-2020 15:55

**Equipment Name** : HIMEDIA BOD SMH-33 673/663

**Temperature (°C)** **Low Alarm** : 0.0 **High Alarm** : 0.0

	Date	Time	Set	-----				Temperature (°C)	-----
Rec. No	dd-MM-yyyy	HH:MM	Temp	Ch-1	Ch-2	Ch-3	Ch-4		
61	09-01-2020	14:05	10.0	9.7	10.1	9.8	10.5		
62	09-01-2020	14:10	10.0	9.8	10.1	9.8	10.5		
63	09-01-2020	14:15	10.0	9.9	10.1	10.0	10.5		
64	09-01-2020	14:20	10.0	10.0	10.2	10.0	10.6		
65	09-01-2020	14:25	10.0	10.1	10.3	10.1	10.7		
66	09-01-2020	14:30	10.0	9.9	10.2	10.0	10.6		
67	09-01-2020	14:35	10.0	9.8	10.1	9.8	10.5		
68	09-01-2020	14:40	10.0	9.8	10.1	9.9	10.5		
69	09-01-2020	14:45	10.0	10.0	10.2	10.0	10.6		
70	09-01-2020	14:50	10.0	10.1	10.2	10.0	10.6		
71	09-01-2020	14:55	10.0	10.1	10.3	10.1	10.7		
72	09-01-2020	15:00	10.0	9.9	10.2	10.0	10.6		
73	09-01-2020	15:05	10.0	9.8	10.1	9.9	10.5		
74	09-01-2020	15:10	10.0	9.9	10.2	10.0	10.6		
75	09-01-2020	15:15	10.0	10.1	10.2	10.1	10.6		
76	09-01-2020	15:20	10.0	10.0	10.3	10.1	10.7		
77	09-01-2020	15:25	10.0	9.8	10.2	10.0	10.5		
78	09-01-2020	15:30	10.0	9.8	10.1	9.9	10.6		
79	09-01-2020	15:35	10.0	9.7	10.1	9.8	10.5		
80	09-01-2020	15:40	10.0	9.9	10.1	9.9	10.5		

---

View Data Log Report**Filter From** : 01-01-2020 00:00 **Filter To** : 09-01-2020 23:59**Records From** : 08-01-2020 16:25 **Records To** : 09-01-2020 15:55**Equipment Name** : HIMEDIA BOD SMH-33 673/663**Temperature (°C)** **Low Alarm** : 0.0 **High Alarm** : 0.0

---

	Date	Time	Set	- - - - -				Temperature (°C)	- - - - -			
Rec. No	dd-MM-yyyy	HH:MM	Temp	Ch-1	Ch-2	Ch-3	Ch-4					
81	09-01-2020	15:45	10.0	9.8	10.0	9.9	10.5					
82	09-01-2020	15:50	10.0	10.0	10.2	10.0	10.6					
83	09-01-2020	15:55	10.0	10.0	10.3	10.1	10.6					