

Comparison and Feature between ICDAS3.0 and Megvid

Sr. No	ICDAS3.0	Megvid
1.	Department level installation	Company level installation
2.	Server base installation	Cloud base installation
3.	Multiple department installation is supported	Not required multiple installation
4.	Only chamber data and alarm collection	Chamber and data logger data and alarm can be collected in Megvid.
5.	Parameter transfer and setting based on chamber and software connectivity	Parameter transfer and setting is independent of equipment
6.	Once department is switched it will be application for all option	Plant/department selection is given on option level
7.	Data collection is done through TCP/web service bases	Data is send to Megvid through FTP protocol
8.	Data review/ approve can only be done within the ICDAS3.0 software	Data Review/ Approve Email link can be shared to the privileged user to review/ approve data
9.	Security Question and Answer system is not available for User authentication	Security Question and Answer system is available for User authentication
10.	LDAP User login system is available	LDAP User login system is not available
11.	Language selection is not available	Language selection is available on user Plant and User level
12.	SMS modem can be configure to send Alarm SMS	Bulk SMS service provider link need to be configure to send Alarm SMS
13.	License is Client PC base	License is on equipment bases
14.	On 1 Client service maximum 16 equipment can be connected	There is no limit to connect equipment and License is based on equipment bases
15.	Support only chamber adding	Support adding chamber and data logger.
16.	Mobile Application not Supported	Mobile Application Supported for android and iphone also
17.	HMI login User Level not Available	HMI Login by user level is available
18.	Supervisor email is handled in department	Supervisor email is handled in separate option in Manage Authorization
19.	Data Review and Approve via Mobile App is not Possible	Data Review and Approve Via Mobile Application is Possible